# Ni-NTA Spin Kit Handbook

Ni-NTA Spin Kit

Ni-NTA Spin Columns

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**Table 1. Compatibility of reagents with Ni-NTA**

Reagent	Effect	Comments
<b>Buffer reagents</b>		
Tris, HEPES, MOPS	Buffers with secondary or tertiary amines may reduce nickel ions	Up to 100 mM can be used, however sodium phosphate or phosphate-citrate buffer is recommended
<b>Chelating reagents</b>		
EDTA, EGTA	Strip nickel ions from resin	Up to 1 mM has been used successfully in some cases, but care must be taken
<b>Sulfhydryl reagents</b>		
$\beta$ -mercaptoethanol	Prevents disulfide cross-linkages. Can reduce nickel ions at higher concentration	Up to 20 mM can be used. Do not store resin under reducing conditions
DTT, DTE	At high concentrations (>1 mM) resin may turn reversibly brown due to nickel reduction. Up to 10 mM has been tested and shown not to compromise purification or increase nickel leaching.	Up to 10 mM DTT has been used successfully. Do not store resin under reducing conditions.
<b>Detergents</b>		
Nonionic detergents (Triton®, Tween NP-40, etc.)	Removes background proteins and nucleic acids	Up to 2% can be used
Cationic detergents		Up to 1% can be used
Nonionic detergents (@-OG, DM, DDM, Cymal 6, Apo12 9, NG and others)	Resolubilization and purification of membrane proteins	Up to 2% can be used
Zwitterionic detergents (LDAO, CHAPS, CHAPSO)	Removal of background proteins and nucleic acids; purification of membrane proteins	Up to 1% can be used
Anionic detergents (SDS, sarkosyl)		Not recommended, 0.3% has been in some cases

## Ni-NTA spin kit

Triton X-114	Removes endotoxins	Up to 2% can be used
<b>Denaturants</b>		
GuHCl	Solubilize proteins	Up to 6 M
Urea		Up to 8 M
<b>Amino acids</b>		
Glycine		Not recommended
Glutamine		Not recommended
Glutamic acid		Up to 100 mM
Arginine		Not recommended
Histidine	Binds to Ni-NTA and competes with histidine residues in the 6xHis tag	Can be used at low concentrations (20 mM) to inhibit nonspecific binding and, at higher concentrations (>100 mM), to elute the 6xHis-tagged protein from the Ni-NTA matrix
<b>Other additives</b>		
NaCl	Prevents ionic interactions	Up to 2 M can be used, at least 300 mM should be used
MgCl <sub>2</sub>		Up to 4 M
LiCl <sub>2</sub>		Up to 50 mM
CaCl <sub>2</sub>		Up to 5 mM
MgSO <sub>4</sub>		Up to 1 M
Glycerol	Prevents hydrophobic interaction between proteins	Up to 50%
Ethanol	Prevents hydrophobic interactions between proteins	Up to 20%
Triethanolamine	Prevents hydrophobic interactions	Up to 100 mM
Sorbitol, betaine, ectoine	Prevents hydrophobic interactions	Up to 500 mM
Dextran sulfate	Prevents hydrophobic interactions	Up to 0.05% (w/v)
Imidazole	Binds to Ni-NTA and competes with histidine residues in the 6xHis tag	Can be used at low concentrations (20 mM) to inhibit non specific binding and, at higher concentrations (>100

Ni-NTA spin kit

		mM), to elute the 6xHis-tagged protein from the Ni-NTA matrix
Sodium bicarbonate		Not recommended
Hemoglobin		Not recommended
Ammonium		Not recommended
Citrate		Up to 60 mM has been used successfully

## Troubleshooting Guide

	Comments and suggestions
<b>Protein does not bind to the Ni-NTA Spin column</b>	
a) 6xHis tag is not present	Check expression construct. Sequence ligation junctions to ensure that the reading frame is correct. Check for possible internal translation starts (N-terminal tag) or premature termination sites (C-terminal tag).
b) 6xHis tag is inaccessible	Purify protein under denaturing conditions. Move tag to the opposite end of the protein.
c) 6xHis tag has been degraded	Check that the 6xHis tag is not associated with a portion of the protein that is processed.
d) Binding conditions incorrect	Check pH of all buffers. Dissociation of urea often causes a shift in pH. The pH values should be checked immediately prior to use Ensure that there are no chelating or reducing agents present and that the concentration of imidazole is not too high
<b>Protein elutes in the wash buffer</b>	
a) Wash stringency is too high	Lower the concentration of imidazole or increase the pH slightly.
b) 6xHis tag is partially hidden	Purify under denaturing conditions.
c) Buffer conditions incorrect	Check pH of denaturing wash buffer.
<b>Protein precipitates during purification</b>	
a) Temperature is too low	Perform purification at room temperature.
b) Protein forms aggregates	Try adding solubilization reagents such as 0.1% Triton® X-100 or Tween®-20, up to 20 mM β-ME, up to 2 M NaCl, or stabilizing cofactors such as Mg <sup>2+</sup> . These may be necessary in all buffers to maintain protein solubility.
<b>Protein does not elute</b>	
Elution conditions are too mild (protein may be in an aggregate or multimer form)	Elute with decreased pH or increased imidazole concentration.
<b>Protein elutes with contaminants</b>	
a) Binding and wash conditions	Include 10–20 mM imidazole in the binding not stringent enough and wash buffers
b) Contaminants are associated with tagged protein	Add-mercaptoethanol to a maximum of 20 mM to reduce disulfide bonds.

Ni-NTA spin kit

	Increase salt and/or detergent concentrations, or add ethanol/glycerol to wash buffer to disrupt nonspecific interactions (see Table 1, page 11).
c) Contaminants are truncated	Check for possible internal translation start sites of the tagged protein (C-terminal tag) or premature termination sites (N-terminal tag). Prevent protein degradation during purification by working at 4°C or by including protease inhibitors

## Appendix A: Preparation of Guanidine-Containing

### Samples for SDS-PAGE

Since the fractions that contain GuHCl will form a precipitate when treated with SDS, they must either be diluted with water (1:6), dialyzed before analysis, or separated from the guanidine hydrochloride by trichloroacetic acid (TCA) precipitation. TCA-precipitation: Bring the volume of the samples up to 100  $\mu$ l with water, add an equal volume of 10% TCA, leave on ice 20 min, spin 15 min at 15,000 x g in a microcentrifuge, wash pellet with 100  $\mu$ l of ice-cold ethanol, dry, and resuspend in 1x SDS-PAGE sample buffer (5x SDS-PAGE sample buffer is 0.225 M Tris-Cl, pH 6.8; 50% glycerol; 5% SDS; 0.05% bromophenol blue; 0.25 M DTT). In case there is still any GuHCl present, samples should be loaded immediately after boiling for 7 min at 95°C.

## Appendix B: Buffer Compositions

### Bacterial media and solutions

LB medium	10 g/liter tryptone; 5 g/liter yeast extract 10 g/liter NaCl
LB agar	LB medium containing 15 g/liter agar
Psi broth	LB medium, 4 mM MgSO <sub>4</sub> ; 10 mM KCl
Kanamycin stock solution	25 mg/ml in H <sub>2</sub> O, sterile filter, store in aliquots at -20°C
Ampicillin stock solution	100 mg/ml in H <sub>2</sub> O, sterile filter, store in aliquots at -20°C
IPTG (1 M)	238 mg/ml in H <sub>2</sub> O, sterile filter, store in aliquots at -20°C